

B+ve[®]



**Improving
The Quality Of
Your Life Through
Better Health.**

PRODUCT **CATALOGUE**

omexmedicaltechnology.co.in

About Us

Since 1994, Omex Group of Industries has evolved from a single line of infusion-set parts into a diversified, ISO, CE, and WHO-GMP certified healthcare manufacturer, firmly established across India's leading healthcare networks. With decades of experience, we specialize in sterile, high-quality medical devices, produced in fully integrated facilities that maintain complete control over quality, safety, and performance—from raw material to final packaging. Our extensive product range covers IV therapy, surgical disposables, wound care, orthopaedic supports, respiratory and anesthesia care—trusted daily by doctors, hospitals, and caregivers nationwide.

B+ve® - The Heartbeat of Omex

- ▶ At the core of Omex Medical Technology is our proudly registered brand, B+Ve — a name that stands for more than positivity. It reflects precision, protection, and performance in every drop, every wrap, and every breath.
- ▶ Widely trusted across India, B+Ve is known for its clinical reliability, sterile safety, and doctor-led design, supporting everything from infusion therapy and surgical care to respiratory and anesthesia needs.
- ▶ Every B+Ve product is built to meet the stringent demands of modern healthcare, delivering consistent quality and confidence in every use.
- ▶ Sterile. Safe. Trusted. Always.

OUR MISSION

To design and deliver reliable, clinically validated medical products that improve patient outcomes and make caregivers' work safer and simpler.

OUR VISION

To be a global benchmark for affordable, innovative, and sustainable healthcare solutions, touching lives across every continuum of care.

OUR VALUES

Quality. Integrity. Service. Sustainability. We deliver safe, trusted products through ethical practices, fast support, and eco-conscious manufacturing.



Mr. Arvind Patel
Founder

OUR COMPANIES

OMEX INDUSTRIES : Specialised in manufacturing high-precision IV set components, serving as the backbone of infusion therapy with consistent quality and reliability.

OMEX MEDICAL TECHNOLOGY: Expert in a wide range of surgical disposable products, delivering sterile, single-use solutions trusted across hospitals, clinics, and critical care settings.

OMEX PHARMA : Engaged in the supply of high-quality IV fluids, ensuring safe and effective solutions for essential hospital and clinical care.

For over two decades, Omex has been driven by a single purpose: To raise the standard of care through uncompromising quality and continuous innovation.

Thank you for choosing Omex — together, we advance healthcare, one sterile pack at a time.

INDEX

B+ve INFUCARE RANGE

IV SET REGULAR NON VENTED	06
IV SET REGULAR NON VENTED - Winged	06
OMEX FLOW NON VENTED-Winged	06
IV SET VENTED	07
IV SET VENTED Y PORT LUER LOCK	07
IV SAFETY ULTRA PREMIUM NON VENTED	07
MICRO INFUSION SET NON VENTED	08
MICRO INFUSION SET VENTED	08
IV INFUSION KIT	08
FLOW REGULATOR	09
FLOW REGULATOR WITH VENTED IV SET	09
MEASURED VOLUME SETS 110/150ML	10
B T SET REGULAR	11
B T SET DOUBLE DRIP	11
B T SET VENTED	11
B T SET WITH LUER LOCK	11

B+ve FLOWLINE RANGE

IV CANNULAS	13
3-WAY STOP COCK	13
3-WAY STOP COCK WITH EXTENSION LINE	13
UMBILICAL CATHETER	14
PRESSURE MONITORING LINE	14
EXTENSION LINE	14
SCALP VEIN SET	14
CT COILED PRESSURE CONNECTION LINE	14
SPIRAL EXTENSION LINE	14

B+ve NEPHRO-URO CARE RANGE

LATEX FOLEYS CATHETER	16
100 % SILICON FOLEYS CATHETER	16
TRANSURETHRAL RESECTION SET	17
NELTON / URETHRAL CATHETER	17
MALE EXTERNAL CATHETER	17
URINE BAG REGULAR	18
URINE BAG STANDARD	18
URINE BAG PREMIUM	18
PAEDIATRIC URINE BAG	18
URO METER	19
HEMODIALYSIS BLOOD TUBING SET	19
HEMODIALYSIS CATHETER KIT (DL-TL)	19
AV FISTULA NEEDLE	19

B+ve SURGIBASE RANGE

YANKAUER SUCTION SET	21
THORACIC DRAINAGE CATHETER	21
THORACIC TROCAR CATHETER	21
ABDOMINAL DRAINAGE KIT	21
CLOSE WOUND SUCTION DRAIN UNIT	22
MINI CLOSE WOUND SUCTION DRAIN UNIT	22
UNDER WATER SEAL DRAINAGE BAG	22
CORRUGATED DRAINAGE SHEET	23
SURGICAL BLADE	23
SKIN GRAFTING BLADE	23
CAUTERY PENCIL	23
SKIN STAPLER	23

B+ve GASTRO CARE RANGE

INFANT FEEDING TUBE	25
RECTAL CATHETER	25
RYLE'S TUBE (NG Tube)	25

B+ve ANESTORESP CARE RANGE

SUCTION CATHETER	27
CLOSED SUCTION CATHETER	27
NASAL CANNULA	28
OXYGEN MASK	28
HI-CONCENTRATION MASK (NRBM)	28
NEBULIZER MASK KIT	28
NEBULIZER WITH T-MOUTH-PCS	29
SUCTION TUBE WITH CONNECTOR	29
VENTURY OXYGEN MASK WITH 7-DIAL	29
OROPHARYNGEAL/GUEDEL AIRWAYS	29
NASOPHARYNGEAL AIRWAYS	30
ENDOTRACHEAL TUBES	30
TRACHEOSTOMY TUBE	30
SPINAL NEEDLE	30
T-PCS COMPONENT	31
CATHETER MOUNT	31
BREATHING FILTER	31
T-OXYGENATOR WITH TUBING	31
PLAIN VENTILATOR CIRCUIT	31
SINGLE WATER TRAP VENTILATOR CIRCUIT	32
DOUBLE WATER TRAP VENTILATOR CIRCUIT	32
CO-AXIAL BAIN'S CIRCUIT	32
VENTILATOR KIT	32
NON INVASIVE VENTILATION MASK	33

INDEX

B+ve ANESTORESP CARE RANGE

3-BALL RESPIRATORY EXERCISER	33
SILICON RESUSCITATOR KIT	33
RE-BREATHING BAGS	33

B+ve WOUNDSECURE RANGE

PLASTER OF PARIS BANDAGE	35
COMBINE DRESSING	35
ADHESIVE TAPE (SPOOL TYPE)	35
ADHESIVE TAPE (BAMBOO TYPE)	35
MICROPORE PAPER TAPE	35
CREPE BANDAGE	36
GAMJEE ROLL	36
SOFT ROLL	36
ELASTIC ADHESIVE BANDAGE	36
ABSORBANT COTTON WOOL	36

B+ve MEDIESSENTIALS RANGE

EXAMINATION GLOVES	38
SURGICAL GLOVES	38
PLASTIC GLOVES	38
UMBILICAL CORD CLAMP	38
MUCUS EXTRACTOR	38
CANNULA FIXATOR	38
URINE / STOOL CULTURE BOTTLE	39
KARMAN TYPE CANNULA	39
HOT WATER BAG	39
PATIENT ID BAND	39
BLOOD LANCET	39
TOURNIQUETE BELT	39
ALCOHOL SWAB	40
SKIN PREP RAZOR	40
URINE POT 2IN 1	40
BED PAN	40
ECG/ULTRASOUND GEL	40
ECG ELECTRODE	40
SAFE DELIVERY KIT	41
HIV PROTECTION KIT	41
UNDER PAD	41
SURGICAL GOWN	42
SURGICAL APRON	42
BED SHEET	42
SHOE COVERS	42
MEDICAL CAPS	42
3-PLY MASK	42

B+ve ORTHOSUPPORT ESSENTIALS

BODY BELT AND BRACES

LUMBAR SACRAL BELT	44
MATERNITY BELT	44
ABDOMINAL BELT	44
RIB BELT	44

CERVICAL AIDS

SOFT CERVICAL COLLAR WITH BUTTON	44
HARD CERVICAL COLLAR	44
PHILADELPHIA NECK COLLAR	44

FRACTURED AIDS

ARM SLING BAGGY/TROPICAL	45
SHOULDER IMMOBILIZER	45
ELASTIC SHOULDER IMMOBILIZER	45
CLAVICLE BRACE	45

WRIST AND FOREARM

WRIST BINDER	46
WRIST BRACE WITH THUMB	46
TENNIS ELBOW SUPPORT	46
WRIST SPLINT	46
THUMB SPICA SPLINT	46
ELBOW SUPPORT	46

FINGER SPLINT ORTHOSIS

FINGER Mallet	47
FINGER COT	47
FROG SPLINT	47

KNEE AND ANKLE SUPPORT

KNEE CAP & WITH OPEN PATELLA	47
ELASTIC KNEE SUPPORT	47
SKIN TRACTION KIT	47
COMPRESSION STOCKING	48
ANKLE SUPPORT	48
KNEE IMMOBILIZER	48
ANKLET WITH BINDER	48
ANKLE TRACTION	48
ANKLE BRACE	48



InfuCare Range

B+ve®

The B+ve InfuCare Range is designed to support intravenous therapy, where fluids or medications are delivered directly into the bloodstream. This method ensures rapid absorption and precise control, making it essential in both routine and critical care settings.

IV SET REGULAR

NON VENTED



IS 404

- Streamlined drip chamber
- Plain latex port and smooth roller clamp
- Flexible tubing keeps flow steady and fuss-free.

IV SET REGULAR

NON VENTED - Winged




IS 404 A

- Easy-grip winged chamber,
- Plain latex port and dependable roller clamp
- Supple tubing delivers even, hassle-free infusion.

OMEX FLOW

NON VENTED-Winged



Bulb latex port



Winged chamber

IS 402

- Clear-pack set pairing a winged drip chamber
- Bulb latex port and confident roller clamp
- Kink-safe tubing ensures consistent worry-free delivery.

IV SET VENTED



Vented chamber



Bulb latex port



IS 406

- Vented set featuring an air-balancing chamber for rigid bottles
- Bulb injection port and confident roller clamp
- Kink-safe tubing to deliver steady, worry-free infusions.

IV SET VENTED

Y PORT LUER LOCK



Y injection port



Luer-lock with Needle




IS 420

- Ribbon-packed upgrade featuring a vented chamber
- Y-port, secure luer-lock and glide-easy roller clamp
- Crystal tubing supports accurate, contamination-free delivery.

SAFETY ULTRA PREMIUM

NON VENTED



Bulb latex port



Safety chamber

IS-419

- Ribbon-packed premium set with a safety chamber
- Bulb Latex injection port and finely tuned roller clamp
- Soft tubing offers quiet, worry-free performance.

MICRO INFUSION SET

NON VENTED



Micro drip
non vented

IS-409

- Paediatric micro-infusion set calibrated to ≈ 60 drops/ml
- For low-volume dosing; sharp spike and fine-bore tubing
- Smooth roller clamp to keep flow steady and fuss-free.

MICRO INFUSION SET

VENTED



Micro drip
vented

IS-405

- Vented paediatric micro-infusion set delivering ≈ 60 drops/ml;
- easy-turn Roller clamp keep flow accurate and worry-free.

IV INFUSION KIT

VENTED IV SET, 22 G IV CANNULA, ALCOHOL SWAB, CANNULA FIXATOR, 5 ML SYRINGE



IVF-422

- All-in-one sterile infusion kit combining a vented IV set, intravenous cannula, alcohol swab, cannula fixator and disposable syringe
- Ready for quick, organised bedside initiation of fluid therapy.

FLOW REGULATOR



Y injection port



Dial Flow

IVF-421

- Precision flow control for IV fluids with smooth and accurate adjustment.
- Designed for consistent performance and reliability during infusion therapy.

FLOW REGULATOR WITH VENTED IV SET



IVF-421 A

- Combines a flow regulator, vented Drip chamber, Y-port, and Luer lock for versatile IV access.
- Ensures controlled, safe fluid administration in clinical settings.



Dial Flow



Vented Chamber



Y-site



Luer-lock



MEASURED VOLUME SET

110 ML



Y injection port



Needle free
luer lock



MV-412 A MV-412 B

- Calibrated pediatric infusion set with a 110ml chamber
- Available in two variants: one with a bulb latex injection port for manual access, and Another with an integrated y-port and luer lock for closed-system medication delivery.
- Designed for safe, accurate, and controlled infusion in neonatal and pediatric care.

MEASURED VOLUME SET

150 ML



Bulb Injection Site



D-Clamp



MV-421 A MV-421 B

- Pediatric infusion set with a 150ml calibrated burette chamber offered in variants with bulb latex injection port or Y-port with Luer lock.
- Ensures accurate volume control and safe medication administration in neonatal and pediatric fluid therapy.

BLOOD TRANSFUSION SET

REGULAR



Drip chamber with fluid filter



BT-410

- Standard design with transparent drip chamber and sharp spike.
- Ensures smooth, uninterrupted blood flow in routine transfusions.

BLOOD TRANSFUSION SET

DOUBLE DRIP



Double drip Chamber

BT-411 A

- Dual chamber improves air separation and visual monitoring of flow.
- Enhances transfusion safety, especially in critical care.

All BT Sets include a 200-250 micron fluid filter to trap clots and impurities, ensuring safe blood transfusion.

BLOOD TRANSFUSION SET

VENTED



Vented drip Chamber



BT-411B

- Vented chamber allows compatibility with both glass and plastic bottles.
- Ideal for versatile transfusion setups in varied clinical environments.

BLOOD TRANSFUSION SET


WITH LUER LOCK



Bulb latex with Luer-lock

BT-411C

- Equipped with a secure Luer lock connector for leak-proof performance.
- Enables safe transfusion connection to Luer-compatible ports.



FlowLine Range

B+ve®

The B+ve FlowLine Range focuses on vascular access and fluid management, enabling smooth and controlled administration or withdrawal of fluids. It plays a vital role in maintaining consistent flow during medical procedures, ensuring patient safety and clinical efficiency.

INTRA-VENOUS CANNULA

I.V. Cannula



TECHNICAL SPECIFICATION

Size	Catheter OD (mm)	Catheter Length (mm)	Flow Rate (Water) (ml/min)
16G	1.7	45	200
18G	1.3	45	90
20G	1.1	32	60
22G	0.9	25	36
24G	0.7	19	23
26G	0.6	19	13



Luer-lock
Flash back
chamber



Injection Port



Sharp Needle

IV-418

- Injection port for bolus drug delivery without additional punctures.
- Integrated non-return valve prevents backflow of blood or fluids.
- Flexible wings allow easy fixation and enhanced patient comfort.
- Color-coded hubs ensure quick and accurate size identification.

3-WAY STOP COCK



WS-414

- Allows fluid administration and withdrawal through multiple lines.
- 360° rotating tap with male & female Luer connectors.
- Leak-proof, pressure-resistant construction for safe use.

3-WAY STOP COCK WITH EXTENSION LINE



Sizes: 10, 25, 50, 100, 150, 200 cm

ENL-415

- Combines a 3-way stop cock with flexible extension tubing for improved accessibility.
- Reduces line clutter and ensures safe drug delivery from a distance.

UMBILICAL CATHETER



Size : 4, 5, 6, 7, 8 FG

UC-506

- › Sterile catheter for neonatal vascular access via umbilical vessels
- › Used in NICU for fluid administration, blood sampling, and pressure measurement

PRESSURE MONITORING LINE



Size : 100, 150, 200, 250 cm

HPM-416

- › Designed for accurate transmission of pressure from catheter to monitor.
- › Kink-resistant and non-compliant tubing for reliable pressure readings.

EXTENSION LINE



Size : 100, 150, 200 cm

ET-417

- › Flexible tubing to extend IV lines for drug delivery.
- › Kink-resistant and compatible with standard IV connectors.

SCALP VEIN SET



Available Sizes



18 G 20 G 21 G 22 G 23 G 25 G 27 G

SV-414

- › Winged IV set designed for short-term venous access
- › Colour coded Soft butterfly wings easy size ID and fixation

CT COILED PRESSURE CONNECTION LINE



Size : 150, 200, 250 cm

CTPC-417

- › Coiled tubing for pressure monitoring in diagnostic imaging setups.
- › Tangle-free, compact design ensures smooth performance.

SPIRAL EXTENSION LINE



Size : 150, 200 cm

SHPM-419

- › Spiral-shaped line offers flexibility and prevents kinks.
- › Ideal for patient mobility and long-duration infusions.

Nephro-Uro Care Range

B+ve®

The B+ve Nephro-Uro Care Range addresses essential needs in urinary drainage and renal support. It is designed to assist in the management of catheterization, urine collection, and dialysis procedures, ensuring hygiene, efficiency, and patient comfort in both hospital and home care

LATEX FOLEYS CATHETER

2-WAY , 3-WAY



2 Way



3 Way



Pressure Resistant Balloon



Optimal eyes

Available Sizes

- 6 FR
- 8 FR
- 10 FR
- 12 FR
- 14 FR
- 16 FR
- 18 FR
- 20 FR
- 22 FR


FC-301

FC-302


- Made from high-quality latex with smooth silicone coating for comfort.
- 2-way and 3-way variants available for standard and continuous bladder irrigation use.
- Symmetrical balloon ensures uniform expansion and minimal trauma.

100 % SILICON FOLEYS CATHETER

2-WAY , 3-WAY



2 Way



3 Way



Pressure Resistant Balloon


Available Sizes

- 6 FR
- 8 FR
- 10 FR
- 12 FR
- 14 FR
- 16 FR
- 18 FR

FC-304

- Non-reactive, biocompatible catheter ideal for long-term catheterization.
- Available in 2-way and 3-way options with high fluid compatibility.
- Transparent body for easy observation of flow and cleanliness.

TRANSURETHRAL RESECTION SET TUR SET



TUR-312

- Designed for continuous bladder irrigation in urology procedures.
- Smooth, kink-resistant tubing with large bore for efficient flow.
- Compatible with standard urological irrigation systems.



Y connector



Clamp



Latex tube



Spike

NELTON / URETHRAL CATHETER



Available Sizes

- 6 FG
- 8 FG
- 10 FG
- 12 FG
- 14 FG
- 16 FG
- 18 FG
- 20 FG
- 22 FG

NC-303 UC- 305

- Soft, rounded tip for painless insertion and minimal urethral trauma.
- Sterile, single-use catheter ideal for intermittent catheterization.

MALE EXTERNAL CATHETER



Available Sizes

EX- Small	20 mm
Small	25 mm
Medium	30 mm
Large	35 mm
EX- Large	40 mm

ME-306

- Non-invasive urinary device for male patients with incontinence.
- Secured using provided adhesive tape for a firm and comfortable fit.

URINE BAG REGULAR

OU-102

URINE BAG STANDARD

OU-103 A

URINE BAG PREMIUM


RIBBON PACK

OU-103 B



- Urine collecting bag used for continuous urine drainage.
- Two-litre capacity with super smooth kink-resistant tubing.
- Fitted with an efficient non-return valve and top outlet for backflow prevention.

PAEDIATRIC URINE BAG



OU-104

- Specially designed for pediatric patients with soft, skin-friendly adhesive.
- Graduated chamber provides accurate urine output measurement in Pediatrics.

URO METER



OU-106

- Integrated urometer for precise hourly urine output monitoring.
- Ideal for ICU, critical care, and post-operative patients.
- Features a needle-free sampling port and clear graduated scale.

HEMODIALYSIS BLOOD TUBING SET



BTS-950

- High-quality, sterile tubing designed for use in hemodialysis.
- Kink-resistant and pressure-tolerant to ensure uninterrupted blood flow.
- Compatible with all standard hemodialysis machines.


HEMODIALYSIS CATHETER KIT (DL-TL)



HDC- 1101

- Sterile catheter set available in double-lumen (DL) and triple-lumen (TL) variants.
- Designed for temporary vascular access during hemodialysis.
- Includes catheter, guidewire, dilator, and all necessary components for safe insertion.

AV FISTULA NEEDLE (SINGLE / TWIN PACK)



Sizes: 16, 17 G

AV-1005

- Sharp, thin-wall needle for reliable access to AV fistulas.
- Available in single and twin pack configurations with color-coded wings.
- Features rotating luer lock and flexible wings for easy fixation.

SurgiBase Range

B+ve®

The B+ve SurgiBase Range supports surgical procedures through reliable, sterile, and easy-to-use tools designed for suction, drainage, incision, and operative care. This range plays a key role in maintaining safety, precision, and hygiene throughout the surgical process.

YANKAUER SUCTION SET



Plain tip



Crown tip



YSC-825

- Used for oral and surgical site suctioning during procedures.
- Rigid, transparent handle with smooth atraumatic tip.
- 2.5-m. kink-resistant PVC tubing with universal connectors.
- Available with standard And crown tip designs.

THORACIC DRAINAGE CATHETER



Length: 40-45 cm

Sizes: 12 to 36 FG

CDC-834

- Used For draining air or fluid from the pleural cavity.
- Flexible PVC catheter with radio-opaque line.
- Rounded, multiple eyes for smooth fluid evacuation.

THORACIC DRAINAGE CATHETER WITH TROCAR



Length: 45 cm

Sizes: 12 to 36 FG

TDC-835

- Includes sharp trocar for quick insertion in emergencies.
- Marked catheter with radio-opaque line and multiple side eyes.

ABDOMINAL DRAINAGE KIT



ADK-836

- Post-operative kit used to drain fluid from the abdominal cavity.
- Catheter length: 45 cm, with multiple perforations.
- 2-liter graduated bag with hanger and top outlet.

CLOSE WOUND SUCTION DRAIN UNIT

**CWS-871**

- For continuous, low-vacuum drainage of post-surgical wound fluid.
- Includes graduated bellows for low vacuum suction.
- Provided with curved needle and matching perforated catheter for efficient fluid collection.
- Supplied with spare perforated catheter and connection tube with clamp for secure, controlled drainage.

Bellow

400 ml
800 ml

Available Sizes

- 10 FG
- 12 FG
- 14 FG
- 16 FG
- 18 FG



Connecting tube with clamp



Graduated bellows



Curved needle with catheter



Spare perforated catheter

MINI CLOSE WOUND SUCTION DRAIN UNIT



Bellow 50 ml

Sizes: 6 FG, 8FG

CWS-871A


- Ideal for minor surgeries such as plastic surgery, ENT, and small orthopedic procedures

UNDER WATER SEAL DRAINAGE BAG (ICD BAG)

**WDB-837**

- For safe pleural drainage while preventing air re-entry.
- One-way valve for air-tight seal, Graduated 1000 ml bag with transparent chamber for fluid monitoring.
- Comes with 150 cm kink-resistant tubing and hanger.

CORRUGATED DRAINAGE SHEET



Length: 250 mm

Width: 25 mm

CD-842

- Used for passive drainage of surgical wounds.
- Soft PVC with corrugated design for passive wound drainage.
- Sterile x-ray opaque line and easy to cut as per requirement.

SURGICAL BLADE



Filament 3: 10, 11, 12, 15

Filament 4: 20, 21, 22, 23, 24

SB-871

- Carbon/Stainless steel blade with sharp precision edge.
- gamma sterilized, individually packed in peel-open foil.

SKIN GRAFTING BLADE



Length: 158 mm

SGB-874

- Stainless steel blade for clean, uniform graft harvesting.
- Compatible with standard grafting handles.

CAUTERY PENCIL



CP-873

- For electrosurgical cutting & coagulation during operations.
- Lightweight pencil with push-button activation. Compatible with monopolar cautery machines standard 3-pin connector.

SKIN STAPLER



ST-893

- Lightweight, ergonomic design for precise staple placement
- used in general surgery, orthopedics, and emergency procedures




GastroCare Range

B+ve®

The B+ve GastroCare Range is designed to support gastrointestinal access and feeding, particularly in patients unable to take nutrition orally. It plays a vital role in procedures like nasogastric feeding, gastric decompression, and rectal administration—ensuring comfort, safety, and clinical efficiency.

INFANT FEEDING TUBE



Available Sizes

- 5 FG
- 6 FG
- 7 FG
- 8 FG
- 9 FG
- 10 FG

Length: 52 cm

IF-501 IF-502

- Used for enteral feeding in neonates and infants.
- Soft, non-toxic, medical-grade PVC for gentle insertion.
- Marked tube with radio-opaque line, closed distal end with two lateral eyes.
- Available in : Pocket pack and straight pack

RECTAL CATHETER




Length: 40 cm

Sizes: 18 to 30 FG

RC-307

- Used for rectal administration of medication or enema.
- Smooth, rounded tip with twin lateral eyes for optimal flow.
- Made from non-irritant, medical-grade PVC.

RYLE'S TUBE (NG Tube)



Length: 105 cm

RT-503

- For gastric decompression, feeding, or medication administration.
- Made from soft Medical gradePVC with X-ray opaque line and graduation marks.
- Four lateral eyes, closed rounded tip, distal metal olive for easy navigation.



Stopper Connector



Closure Connector

Available Sizes

- 8 FG
- 10 FG
- 12 FG
- 14 FG
- 16 FG
- 18 FG
- 20 FG




AnestoResp Care Range

B+ve®

The B+ve AnestoResp Care Range combines respiratory therapy and anesthesia support, focusing on airway management, oxygen delivery, and ventilation. It is essential during surgical procedures and critical care, helping maintain patient breathing, sedation, and respiratory function with safety and precision.

SUCTION CATHETER



Length: 53 cm

SC-601 SC-607 SC-608

- Soft, non-toxic catheter for airway suctioning.
- Rounded open tip with two lateral eyes; color-coded hub.
- Thumb control valve for precise suction regulation.
- Available in pocket & straight pouch



Thumb control connector



Plain connector

Available Sizes

- 6 FG
- 8 FG
- 10 FG
- 12 FG
- 14 FG
- 16 FG
- 18 FG
- 20 FG

CLOSED SUCTION CATHETER



Polymer sheath



Suction knobs



Washing port with Shut of valve

Available Sizes

S.No.	Size	OD of Catheter	ID of Catheter	Effective Length
1	FG-6	2.00 mm	1.40 mm	300 mm
2	FG-8	2.70 mm	1.80 mm	300 mm
3	FG-10	3.30 mm	2.40 mm	600 mm
4	FG-12	4.00 mm	2.80 mm	600 mm
5	FG-14	4.70 mm	3.20 mm	600 mm
6	FG-16	5.30 mm	3.80 mm	600 mm

CSC-1102

- Allows tracheal suctioning without disconnecting the ventilator, minimizing infection risk.
- Comes with swivel elbow, T-piece, and optional MDI port for bronchodilator delivery.
- Soft, transparent catheter with depth markings and protective sleeve to prevent contamination.
- with irrigation port for catheter flushing.

NASAL CANNULA




Sizes: Adult, Neonatal & Pediatric

NC-818

- Dual-pronged nasal interface for low-flow oxygen delivery.
- Soft PVC tubing: 2 m, kink-resistant with universal connector.

OXYGEN MASK



Sizes: Adult & Pediatric

OXY-605

- Transparent, anatomical mask with adjustable elastic strap.
- Soft PVC mask with aluminum nose clip.
- Tubing: 2 m length, Sizes: Adult & Pediatric.

HI-CONCENTRATION MASK (NRBM)



Sizes: Adult & Pediatric

HCM-817

- Non-rebreathing mask with reservoir bag.
- Provides up to 99% FiO₂, one-way valves for safety.
- Reservoir bag : 600–1000 ml.

NEBULIZER MASK KIT



Sizes: Adult & Pediatric

NB-604

- Designed for aerosol therapy to deliver medication directly to the lungs.
- Includes 6–8 ml nebulizer chamber and 2 m oxygen tubing for optimal mist delivery.

NEBULIZER WITH T-MOUTH-PCS



NB-606

- T-piece with mouthpiece, nebulizer jar & oxygen tubing.
- Delivers medication without using a face mask.
- Chamber capacity: 6–8 ml, tubing: 2 m.

SUCTION TUBE WITH CONNECTOR



ST-826

- Kink-resistant tubing with universal connectors.
- For OR and ICU suction procedures.

VENTURY OXYGEN MASK WITH 7-DIAL



Sizes: Adult

HCM-817

- Delivers fixed oxygen concentration (24% to 60%).
- 7-dial color-coded adapter with clearly marked % levels.

OROPHARYNGEAL/ GUEDEL AIRWAYS



Available Sizes



GA-816

- Oropharyngeal airway with color-coded bite block.
- Smooth, curved shape for easy insertion.

NASOPHARYNGEAL AIRWAYS



Sizes: 12-36 FR

NA-868

- Soft, flexible airway for nasal insertion.
- Flared end prevents displacement.

ENDOTRACHEAL TUBES CUFFED/UNCUFFED



Sizes: 2.0 mm to 8.5 mm ID

ET-602

ET-603

- High-volume, low-pressure cuff (for cuffed type).
- Murphy eye & radio-opaque line; standard 15 mm connector.

TRACHEOSTOMY TUBE CUFFED/UNCUFFED



Sizes: 4.0 to 8.0 mm ID.

TCS-607

- Curved tube with standard 15 mm connector.
- Available in cuffed and plain variants.

SPINAL NEEDLE



Length 90 mm

SPN-609

- Features a Quincke-type cutting bevel for smooth dura penetration.
- Includes a stylet to Prevent needle blockage during insertion.

Available Sizes

- 18 G
- 20 G
- 21 G
- 22 G
- 23 G
- 24 G
- 25 G
- 26 G
- 27 G

T-PCS COMPONENT



- › Multi-access connector used in respiratory and nebulization therapy
- › Allows integration of nebulizer chamber with oxygen supply or ventilator circuit
- › Ensures uninterrupted airflow during medication delivery

CATHETER MOUNT



CM-819

- › Flexible link between airway device and breathing circuit.
- › With double swivel elbow and Luer port for EtCO₂/suction access.
- › Universal 15 mm/22 mm connectors, available in corrugated type.

BREATHING FILTER (BVF/HME)



BF-826

- › Filters for preserving airway humidity and filtering pathogens
- › Reduces risk of cross-contamination in ventilated patients

T-OXYGENATOR WITH TUBING (T- Recovery kit)



TO-820

- › Used for oxygen delivery during recovery or post-extubation.
- › T-shaped connector with exhalation port and 2 m oxygen tubing.
- › Ideal for tracheostomy patients or spontaneous breathers requiring humidified oxygen.

PLAIN VENTILATOR CIRCUIT



Sizes: Adult, Pediatric

VC-822

- › Used for general mechanical ventilation in ICU, OT, and emergency care.
- › Compatible with all standard ventilator machines.
- › Made from corrugated tubing with universal 22 mm connectors.

SINGLE WATER TRAP VENTILATOR CIRCUIT




Sizes: Adult, Pediatric

VC-823

- Includes a water trap to collect condensate and maintain airflow.
- Suitable for short to moderate-duration ventilation support.

DOUBLE WATER TRAP VENTILATOR CIRCUIT



Sizes: Adult, Pediatric

VC-824

- Includes two water traps for enhanced condensate management and infection control.
- Ideal for long-term mechanical ventilation in critical care.

CO-AXIAL BAIN'S CIRCUIT



BHV-867

- A co-axial breathing system used for delivering anesthesia gases.
- Includes 22 mm outer and 6 mm inner tubing, Length: 1.6 m, with APL valve and rebreathing bag.

VENTILATOR KIT

(VENTILATOR CIRCUIT WITH HME FILTER AND CATHETER MOUNT)



Sizes: Adult, Pediatric

VK-870

- Ready-to-use kit for mechanical ventilation, includes plain ventilator circuit, HME filter, and catheter mount.
- Provides secure airway connection, filtration, and humidification in ICU and anesthesia care.

NON INVASIVE VENTILATION MASK



Size : Small, Medium, Large

NIV-869

- Full-face NIV mask for use with BiPAP/CPAP systems.
- Vented (with exhalation ports) and non-vented types.
- Made from soft, medical-grade silicone for a secure seal.
- Standard 22 mm connector with adjustable headgear.

3-BALL RESPIRATORY EXERCISER



BS-828

- Visual incentive spirometer with 3 color-coded balls.
- Encourages deep breathing to improve lung function.
- Ball resistance levels: 600, 900, 1200 cc/sec.

SILICON RESUSCITATOR KIT (AMBU BAGS)



Sizes: Adult, Pediatric, Neonatal

AB-829

AB-830

AB-831

- Manual resuscitation device used for emergency respiratory support.
- Includes self-inflating bag, oxygen reservoir bag, oxygen tubing, and anatomically shaped mask.
- Equipped with pressure limiting valve (in pediatric models) to avoid lung overinflation.

RE-BREATHING BAGS



Capacity : 0.5, 1, 2 Ltr.

RB-840

- Reservoir bag for anesthesia breathing circuits.
- 22 mm neck for secure attachment.



WoundSecure Range

B+ve®

The B+ve WoundSecure Range is designed for effective wound care and post-surgical recovery. It supports protection, absorption, and healing by maintaining a clean and sterile environment—helping reduce the risk of infection and promoting faster tissue repair.

PLASTER OF PARIS BANDAGE



Width : 7.5 cm, 10 cm, 15 cm

Length: 2.7 m, 3 m

PPB-900 PPB-901 PPB-902

- Used for orthopedic immobilization in fractures and sprains.
- Fast-setting POP on gauze roll; provide from support after application.

COMBINE DRESSING



Sizes: 10x10 cm, 10x20 cm

CD-905 CD-906

- Sterile absorbent dressing for post-operative and heavy wound drainage.
- Non-woven outer with high-absorbency cotton filling.

ADHESIVE TAPE (SPOOL TYPE)



Width : 2.5 cm, 5 cm, 7.5 cm, 10 cm

Length: 5 m

AT-909

- Zinc oxide adhesive tape on plastic spool.
- Used for securing dressings and catheters.

ADHESIVE TAPE (BAMBOO TYPE)



Width : 2.5 cm, 5 cm, 7.5 cm, 10 cm

Length: 5 m, 8 m, 9 m, 10 m

AT-910

- Strong, breathable adhesive tape with a bamboo type pack
- Provides firm fixation of dressings and medical devices

MICROPORE PAPER TAPE



Width : 1.25 cm, 2.5 cm, 5 cm, 7.5cm

Length: 5 m and 9.14 m

PT-972

- Hypoallergenic, gentle tape ideal for sensitive skin.
- Breathable and easy to tear by hand.

CREPE BANDAGE



Width : 6 cm, 8 cm, 10 cm, 15 cm

Stretched Length: 4 m

CCB-913 - 916

- Elastic support bandage for compression and muscle support.
- Washable, re-usable, and breathable.

GAMJEE ROLL



Width : 10 cm, 15 cm

Length: 3 m

GR-911 GR-912

- Highly absorbent cotton padding with gauze covering.
- Used for dressing wounds and absorbing exudate.

SOFT ROLL/COTTON CAST PADDING



Width : 10 cm, 15 cm

Length: 3 m

SR-930 SR-931

- Non-woven padding roll for use under plaster or synthetic casts.
- Prevents pressure sores and provides cushioning.

ELASTIC ADHESIVE BANDAGE



Width : 6 cm, 8 cm, 10 cm, 15 cm

Stretched Length: 4/6 m

EB-919 - 925

- Stretchable adhesive bandage for compression, support, and dressing fixation.
- Strong adhesion with breathable fabric.

ABSORBANT COTTON WOOL



Net / Gross weight 50, 100, 200, 400, 500 gm.

ACW-957-961

- Used for general wound cleaning, padding, and absorption.
- Made from pure, bleached cotton with high absorbency.

MediEssentials Range

B+ve®

The B+ve MediEssentials Range includes everyday clinical support items used across wards, clinics, and diagnostic settings. These essentials play a vital role in maintaining hygiene, improving workflow, and ensuring readiness for a wide range of medical procedures.

EXAMINATION GLOVES



Size : Small, Medium, Large, X-Large.

EGL-802

- › Ambidextrous, non-sterile, gloves for routine examination and general use.
- › Made from latex / nitrile, powder-free or powdered.

SURGICAL GLOVES



Size : 6.0, 6.5, 7.0, 7.5, 8.0.

SGL-801

- › Sterile paired gloves for surgical procedures with anatomical left-right fit.
- › Made from latex, powdered with micro-roughened surface for better grip.

PLASTIC GLOVES



PEG-803

- › Transparent, loose-fit gloves for food handling, cleaning, and medical use.
- › Made from LDPE, single-use, economical option.


UMBILICAL CORD CLAMP



UCC-501

- › Used for clamping and cutting the newborn's umbilical cord.
- › Features a secure locking system to prevent accidental reopening.

MUCUS EXTRACTOR



ME-601

- › For aspiration of mucus in neonates to ensure clear airway.
- › Transparent chamber with suction tube and mucus trap.

CANNULA FIXATOR



ICF-927

- › Securely holds IV cannula in place to prevent displacement.
- › Hypoallergenic adhesive base with soft foam padding.

URINE / STOOL CULTURE BOTTLE



Size : 30, 60, 100, 120ml

CB-832

- › Sterile containers for specimen collection and microbial culture.
- › Leak-proof, tamper-evident cap with labelling space.

KARMAN TYPE CANNULA



Size : 4 mm to 12 mm

KC-835

- › Used for uterine aspiration or evacuation.
- › Soft, flexible PVC cannula with rounded, open end.
- › Compatible with suction apparatus.


HOT WATER BAG



HW-885

- › Reusable hot water bag for heat therapy and pain relief
- › with secure screw cap
- › Ideal for cramps, muscle aches, and joint stiffness

PATIENT ID BAND



PIB-886

- › for patient tagging and tracking
- › easy to write or print on
- › Adjustable size for comfortable, secure fit

BLOOD LANCET



BL-873

- › Sterile, single-use with ultra-sharp tip for painless puncture.

TOURNIQUETE BELT



TQB-892

- › Used to temporarily restrict blood flow during venipuncture or IV insertion
- › Quick-release mechanism ensures easy application and removal

ALCOHOL SWAB



AS-863

- Pre-saturated with 70% isopropyl alcohol.
- Used to disinfect skin prior to lancing or injection.

SKIN PREP RAZOR



SB-859

- Single-use razor for pre-surgical hair removal.
- Fixed stainless-steel blade with protective design.

URINE POT 2 IN 1



OU-107

- Design for male and female use
- Made from durable, easy-to-clean medical-grade plastic
- Ideal for bedridden or mobility-restricted patients.

BED PAN



OU-108

- Used for toileting of bedridden or immobile patients
- Ergonomic design ensures comfort and ease of use

ECG/ULTRASOUND GEL

250ML/5 LTR



BL-873

- Conductive gel for ECG, ultrasound, and physiotherapy procedures.
- Non-greasy, water-soluble, hypoallergenic formula.

ECG ELECTRODE



ECGE-892

- Pre-gelled, self-adhesive electrode with high-conductivity hydrogel and snap connector for reliable ECG signal acquisition

SAFE DELIVERY KIT



DK-887

- PVC Drape, PVC Apron/Gown, Disposable Delivery Sheet with hole, Cord Clamp, Hole towel for episiotomy suturing, Wrapper for a newborn baby, Mask, Cap, Shoe Cover, Gloves
- Complete sterile kit for hygienic childbirth.

HIV PROTECTION KIT



HK-848

- PPE kit for safe handling in HIV/suspected exposure cases.
- Includes: Gown, Mask, Cap, Shoe Cover, Gloves, Goggles, Waste Bag.

UNDER PAD



Sizes: 60×60 cm and 60×90 cm

OU-107

- Disposable underpad with absorbent core and leak-proof backing
- Ideal for incontinence care and surface protection in medical settings

SURGICAL GOWN**AS-863**

- › Sterile, fluid-repellent SMS gown with knitted cuffs and waist ties for full-body barrier protection.
- › available in Various Sizes

SURGICAL APRON**SB-859**

- › Disposable, full-length apron available in breathable non-woven or waterproof LDPE variants, offering fluid-resistant front-body protection with tie-back closure

BED SHEET**DK-887**

- › Disposable, hypoallergenic bed sheet made from soft spun-bond non-woven fabric for hygienic patient protection; fluid-repellent, lint-free

SHOE COVERS**DSC-814**

- › Made of LDPE or non-woven material, elastic-fit, anti-skid.
- › Universal size, one-time use.

MEDICAL CAPS**DC-805**

- › Non-woven, breathable, with elastic stretch fit.

3-PLY MASK**DFM-810**

- › Includes Melt-blown filter layer, adjustable nose pin, and soft ear loops.

OrthoSupport Essentials

B+ve®

The B+ve OrthoSupport Essentials range provides targeted support for musculoskeletal injuries, post-surgical recovery, and joint stabilization. Designed for comfort and functionality, it helps reduce pain, restrict movement, and promote healing in orthopedic care.

BODY BELTS AND BRACES

LUMBAR SACRAL BELT



Size : S, M, L, XL, XXL

LP-948

- Provides firm lower back support and compression

MATERNITY BELT



Size : Universal

MB-964

- Offers abdominal and back support during pregnancy

ABDOMINAL BELT



Size : S, M, L, XL, XXL

AB-941

- Supports post-surgical abdominal muscles

RIB BELT



Size : S, M, L, XL, XXL, Uni., Child

RB-926

- Provides gentle circumferential support for rib fractures and chest injuries

CERVICAL AIDS

PHILADELPHIA NECK COLLAR



Size : S, M, L, XL

PPC-968

- Rigid, two-piece cervical collar for immobilization of the neck in trauma or post-surgical care

SOFT CERVICAL COLLAR WITH BUTTON



Size : S, M, L, XL

CC-944

- Provides mild neck support and warmth for strains and sprains; breathable foam construction

HARD CERVICAL COLLAR



Size : S, M, L, XL

CC-945

- Offers rigid cervical immobilization after injury or surgery; two-piece plastic design with tracheal opening

FRACTURED AIDS

ARM SLING
BAGGY/TROPICAL

Size : S, M, L, XL, XXL, Child

AS-928

- › Supports the forearm and shoulder after fractures or sprains
- › adjustable strap

SHOULDER IMMOBILIZER



Size : Universal

SI-929 A

- › Restricts shoulder movement for post-injury stabilization

ELASTIC SHOULDER
IMMOBILIZER

Size : S, M, L, XL, XXL

SI-929 B

- › Elastic wrap-around design ensures secure and adjustable compression

CLAVICLE BRACE



Size : S, M, L, XL, XXL, Child

CB-931

- › Figure-of-eight brace to align and support fractured clavicles;

WRIST AND FOREARM

WRIST BINDER



Size : Universal

WS-935

- › Elastic wrap that delivers mild, adjustable compression for general wrist strains

WRIST BRACE WITH THUMB



Size : Universal

WG-934

- › Rigid-stay brace stabilizing wrist and thumb for sprains or De Quervain's

TENNIS ELBOW SUPPORT



Size : Universal

TE-938

- › Adjustable strap with pressure pad to relieve lateral epicondylitis pain

WRIST SPLINT



Size : Universal

WAS-936

- › Supports and immobilizes the wrist in cases of sprain, strain, or post-cast recovery
- › Adjustable straps ensure a snug, comfortable fit

THUMB SPICA SPLINT



Size : Universal

KS-972

- › Immobilizes the thumb while allowing free wrist and finger movement
- › Provides support in cases of thumb injuries, ligament strain, or post-fracture recovery

ELBOW SUPPORT



Size : S, M, L

STK-955

- › Provides compression and support for elbow strains, sprains, or arthritis
- › Helps reduce pain and inflammation

FINGER SPLIT ORTHOSIS

FINGER MALLET



Size : S, M, L

MFS-939

- Rigid plastic splint that immobilizes the DIP joint in mallet injuries

FINGER COT



Size : S, M, L

FC-958

- Soft protective sleeve shielding fingertip wounds from moisture and impact

FROG SPLINT



Size : S, M, L

FS-959

- Aluminium-reinforced "frog" design to stabilize fractured or sprained fingers

KNEE AND ANKLE SUPPORT

KNEE CAP & WITH OPEN PATELLA



Size : S, M, L, XL, XXL

KC-951 KC-952

- Elastic tubular sleeve delivering uniform compression for mild knee strains

ELASTIC KNEE SUPPORT



Size : S, M, L, XL, XXL

KS-972

- Wrap-around elastic brace offering adjustable compression and side stability

SKIN TRACTION KIT



Size : Universal

STK-955

- Complete kit with foam straps, spreader and weight bag for temporary lower-limb traction

KNEE AND ANKLE SUPPORT

COMPRESSION STOCKING

(ABOVE / BELOW KNEE)



Size : S, M, L, XL

CS-975

- Graduated compression stocking to boost venous return and DVT prevention

ANKLE SUPPORT



Size : S, M, L, XL

AS-954

- Below-knee version offering the same therapeutic compression

KNEE IMMOBILIZER



Size : S, M, L, XL, XXL, 3XL

KI-971

- Three-panel splint that locks the knee in full extension after injury or surgery

ANKLET WITH BINDER



Size : Universal

AB-940

- Elastic figure-8 wrap delivering adjustable compression for sprains

ANKLE TRACTION



Size : Universal

AT-932

- Foam stirrup set with cord and weight pouch for lower-limb traction therapy

ANKLE BRACE



Size : Universal

AB-963

- Lace-up brace with side stays providing firm medial-lateral ankle support



B+ve[®]



**Delivering quality
through excellence**

**Delivering quality
through excellence**



Omex Medical Technology

At : Kotadi, Plot No. 1, Shivalik Industrial Park,
Mehsana Highway Ta : Vijapur, Gujarat (India)

✉ sales@omexmedical.in

🌐 www.omexmedicaltechnology.co.in

📷 OMEXMEDICALTECHNOLOGY

📘 OMEX MEDICAL TECHNOLOGY